



Using LOINC codes to build MIDRC cohorts

Paul Kinahan
Vice-Chair, Radiology Research
University of Washington

Funding from NIH-NIBIB Contract 75N92020D00021

Example use case

- You want a cohort of chest CT images, without contrast agent, to develop or test a model that detects, segments, and analyzes nodules to, say, classify a disease or make a prediction.
- You start with the MIDRC data explorer data.midrc.org/explorer
- Under the “Imaging studies / Study Properties” tab you select “CT” for Study Modality and “CHEST” for Body Part Examined in the Study Property filter
- You are (currently) left with 4,300 Imaging studies (aka “Exams”)
- Then you look under Study Description to rule out those with contrast and you get...

[illegible]

Study Description

☐ no data

☐ CT CHEST ANGIOGRAM AND
PULMONARY ARTERIES WITH IV
CONTRAST

☐ CT CHEST W CONTRAST

☐ CT CHEST WO CONTRAST

☐ CT CHEST WITHOUT IV CONTRAST

☐ Thorax^GH_ThoraxRoutine (Adult) 53

☐ CTA CHEST PULMONARY
ANGIOGRAPHY W CONTRAST 52

☐ CT THORAX PE EXAM 49

☐ CT CHEST WO IV CONT 43

☐ CT CHEST W IV CONTRAST 36

☐ CT Thorax W/Contrast 22

☐ CT ANGIOGRAM PULMONARY W
CONTRAST 13

☐ CTA CHEST PULMONARY EMBOLISM 13

☐ CT PE CHEST 12

☐ CT CHEST ANGIOGRAM WITH IV
CONTRAST 11

☐ CT CHEST ANGIO W OR W AND WO IV 10



CT FOR PE 9

ROTOCOL 9

PE (Adult) 9

t) 8

830

MONARY ANGIOGRAM 8

CONTRAST (PETCT) 8

CHEST_WO (Adult) 8

700

ST_WO (Adult) 8

OMEN PELVIS WO 7

IV CONTRAST LUNG 7

265

ontrast 7

7

ST_WITHOUT (Adult) 7

265

Adult) 7

LVIS WITH IV CONT 6

141

CT CHEST 6

CT CHEST (S) 6

☐ CT CHEST ANGIOGRAPH W IV
CONTRAST PULMONAR EMBOLISM 6

☐ CT CHEST LUNG CANCER SCREENING
WITHOUT CONTRAST ANNUAL 6

☐ CT CHEST WO IVCON 6

☐ JCCI CT CHEST ENHANCED 6

☐ THORAX/ABDOMEN/PELVIS + CONT 6

☐ CT CORONARY CALCIUM SCORE 3

☐ CT CT PULMONARY ANGIOGRAPHY 3

☐ CT Chest/Abd/Pelvis W/Contrast 3

☐ CT LUNG NODULE F/U 3

☐ CTA CHEST FOR PE 3

☐ JHN CHEST W/ OP 3

☐ Thorax^CT_CHEST_WITH (Adult) 3

☐ Thorax^DE_PE (Adult) 3

☐ Thorax^FLASH_CAP_WITH (Adult) 3

☐ ANGIO - CHEST (PE STUDY) 2

☐ CAP OP 2

☐ CHEST IV IP 2

☐ CHEST WITH(Adult) 2

☐ CT ABDOMEN PELVIS W CONTRAST 2

☐ CT ANGIOGRAM CHEST PULMONARY
EMBOLISM W CONTRAST 2

☐ CT ANGIOGRAM PULMONARY
ABDOMEN AND PELVIS W CONTRAST 2

☐ CT ANGIOGRAPHY CHEST 2

☐ CT CARD W/CON PULMNARY VEIN 2

☐ CT CARDIAC CHEST W/ (ROBOTIC &
CORONARY) 2

☐ CT CHEST ABDOMEN PELVIS WO
CONTRAST (ROUTINE) 2

☐ CT CHEST PULMONARY ANGIOGRAM 2

☐ CT CHEST W CON 2

Study Description



- | | | |
|--------------------------|--|-----|
| <input type="checkbox"/> | no data | 830 |
| <input type="checkbox"/> | CT CHEST ANGIOGRAM AND PULMONARY ARTERIES WITH IV CONTRAST | 700 |
| <input type="checkbox"/> | CT CHEST W CONTRAST | 265 |
| <input type="checkbox"/> | CT CHEST WO CONTRAST | 265 |
| <input type="checkbox"/> | CT CHEST WITHOUT IV CONTRAST | 141 |
| <input type="checkbox"/> | CTA PE CHEST W | 111 |
| <input type="checkbox"/> | CT CHEST PULMONARY EMBOLISM W CONTRAST | 99 |
| <input type="checkbox"/> | 01_PM_Thorax_Plain(Adult) | 98 |
| <input type="checkbox"/> | Thorax^CT_PE (Adult) | 74 |
| <input type="checkbox"/> | CT CHEST WITH CONTRAST | 69 |
| <input type="checkbox"/> | CT CHEST WO | 68 |
| <input type="checkbox"/> | CT CHEST WITH IV CONTRAST | 63 |
| <input type="checkbox"/> | CT CHEST WITHOUT CONTRAST | 59 |
| <input type="checkbox"/> | CT CHEST WO IV CONTRAST | 59 |
| <input type="checkbox"/> | CT CHEST W/O CONTRAST | 58 |
| <input type="checkbox"/> | CT CHEST PULMONARY EMBOLISM (CTPE) | 55 |
| <input type="checkbox"/> | Thorax^GH_ThoraxRoutine (Adult) | 53 |
| <input type="checkbox"/> | CTA CHEST PULMONARY ANGIOGRAPHY W CONTRAST | 52 |
| <input type="checkbox"/> | CT THORAX PE EXAM | 49 |
| <input type="checkbox"/> | CT CHEST WO IV CONT | 43 |
| <input type="checkbox"/> | CT CHEST W IV CONTRAST | 36 |
| <input type="checkbox"/> | CT Thorax W/ Contrast | 30 |

- | | | | | | |
|--------------------------|--------------------------------------|----|--------------------------|----------------------------------|----|
| <input type="checkbox"/> | CTA CHEST (PE STUDY) W CONTRAST | 27 | <input type="checkbox"/> | CT CHEST WITHOUT CONTRAST (S) | 10 |
| <input type="checkbox"/> | CTA CHEST PULMONARY EMBOLUS W/ IVCON | 24 | <input type="checkbox"/> | CT ANGIO CHEST FOR PE | 9 |
| <input type="checkbox"/> | CT Angio Pulmonary | 23 | <input type="checkbox"/> | CT CHEST PE PROTOCOL | 9 |
| <input type="checkbox"/> | CT CHEST HIGH RESOLUTION | 23 | <input type="checkbox"/> | Thorax^FLASH_PE (Adult) | 9 |
| <input type="checkbox"/> | CT CHEST ABDOMEN PELVIS W CONTRAST | 22 | <input type="checkbox"/> | CHEST WO(Adult) | 8 |
| <input type="checkbox"/> | CT | | <input type="checkbox"/> | CT CHEST PULMONARY ANGIOGRAM | |
| <input type="checkbox"/> | CT/ | | <input type="checkbox"/> | Thorax^DE_PE (Adult) | |
| <input type="checkbox"/> | JCC | | <input type="checkbox"/> | Thorax^FLASH_CAP_WITH (Adult) | |
| <input type="checkbox"/> | CT W/ | | <input type="checkbox"/> | ANGIO - CHEST (PE STUDY) | |
| <input type="checkbox"/> | THC | | <input type="checkbox"/> | CAP OP | |
| <input type="checkbox"/> | Thc | | <input type="checkbox"/> | CHEST IV IP | |
| <input type="checkbox"/> | CT | | <input type="checkbox"/> | CHEST WITH(Adult) | |
| <input type="checkbox"/> | CT IV (| | <input type="checkbox"/> | CT ABDOMEN PELVIS W CONTRAST | |
| <input type="checkbox"/> | CT | | <input type="checkbox"/> | CT CHEST WO IVCON | 6 |
| <input type="checkbox"/> | CT | | <input type="checkbox"/> | JCCI CT CHEST ENHANCED | 6 |
| <input type="checkbox"/> | Thc (Ad | | <input type="checkbox"/> | THORAX/ABDOMEN/PELVIS + CONTRAST | 6 |
| <input type="checkbox"/> | CT CO | | | | |
| <input type="checkbox"/> | CT/ | | | | |
| <input type="checkbox"/> | CT | | | | |
| <input type="checkbox"/> | CT CHEST ANGIOGRAM WITH IV CONTRAST | 11 | | | |
| <input type="checkbox"/> | CT CHEST ANGIO W OR W AND WO IV | 10 | | | |

- | Procedure | Frequency |
|--|-----------|
| <input type="checkbox"/> CT CHEST, ABDOMEN & PELVIS W CONTRAST | 3 |
| <input type="checkbox"/> CT CORONARY CALCIUM SCORE | 3 |
| <input type="checkbox"/> CT CT PULMONARY ANGIOGRAPHY | 3 |
| <input type="checkbox"/> CT Chest/Abd/Pelvis W/Contrast | 3 |
| <input type="checkbox"/> CT LUNG NODULE F/U | 3 |
| 3 | 3 |
| WITH (Adult) | 3 |
| 3 | 3 |
| WITH (Adult) | 3 |
| 2 | 2 |
| STUDY) | 2 |
| 2 | 2 |
| ST W CONTRAST | 2 |
| 2 | 2 |
| EST PULMONARY
FAST | 2 |
| MONARY
IS W CONTRAST | 2 |
| 2 | 2 |
| CHEST | 2 |
| LMNARY VEIN | 2 |
| 2 | 2 |
| W/ (ROBOTIC & | 2 |
| <input type="checkbox"/> CT CHEST ABDOMEN PELVIS WO CONTRAST (ROUTINE) | 2 |
| <input type="checkbox"/> CT CHEST PULMONARY ANGIOGRAM | 2 |
| <input type="checkbox"/> CT CHEST W CON | 2 |

What is with all the study descriptions?

- There are (in this example) over 270 different study descriptions, some are obscure
- These arise as every hospital uses what it wants, and in general this works okay for local use, but not for external users of pooled data
- Is there a standard?
 - Not one that is widely used
 - However. LOINC (Logical Observation Identifiers Names and Codes) does provide a standard for laboratory test orders and results
 - MIDRC is using the Radlex Playbook / LOINC subset of codes
 - LOINC codes are unique with an algorithmically derived Long Common Name that is designed to look like a Study Description (for this subset)

LOINC Structure

- LOINC Code:
 - Unique numerical identifier
- Long Common Name:
 - Also unique, human readable
- Either field may be used

LOINC

SearchLOINC

LOINC

ct chest class:rad system:chest

Search

RESULTS26

DISPLAYING2001-26

FILTER

VIEWListCard

EXPORT

Status	LOINC	Long Common Name	Component	Property	Timing	System	Scale	Method	Class	Type	Example UCUM Units	Order/Observation
	24627-2	CT Chest	Multisection	Find	Pt	Chest	Doc	CT	RAD	U		Both
	36089-1	CT Chest limited	Multisection limited	Find	Pt	Chest	Doc	CT	RAD	U		Both
	87279-6	CT Chest for screening	Multisection for screening	Find	Pt	Chest	Doc	CT	RAD	U		Both
	29252-4	CT Chest WO contrast	Multisection^WO contrast	Find	Pt	Chest	Doc	CT	RAD	U		Both
	35895-2	CT Guidance for biopsy of Chest	Guidance for percutaneous biopsy	Find	Pt	Chest	Doc	CT	RAD	U		Both
	47366-0	CT Chest limited WO contrast	Multisection limited^WO contrast	Find	Pt	Chest	Doc	CT	RAD	U		Both
	24628-0	CT Chest W contrast IV	Multisection^W contrast IV	Find	Pt	Chest	Doc	CT	RAD	U		Both
	89860-1	CT Chest W inspiration and expiration	Multisection^W inspiration + W expiration	Find	Pt	Chest	Doc	CT	RAD	U		Both

Harmonization of MIDRC Study Descriptions

Input study descriptions

Long-tail problem

Modality	Study Description	Count	Cumulative %
DX	XR CHEST 1 VIEW AP	15945	10.9%
DX	XR CHEST 1 VW, FRONTAL	14137	20.5%
CR	XR PORT CHEST 1V	11941	28.6%
DX	XR CHEST PORTABLE 1 VIEW	9978	35.4%
CR	CHEST PORT 1 VIEW (RAD)-CS	7326	40.4%
CR	XR CHEST 2 VIEWS	6212	44.6%
CR	XR CHEST 1 VW PORTABLE	5873	48.6%
CR	XR CHEST 1 VIEW AP	5335	52.2%
...			
CT	CT P CHEST WO	4	99%
...			
DX	THORACIC SPINE 3 VIEWS	1	100%

Row 444

Row 1385

Harmonized study descriptions

↓ 95%+

Reduction

LOINC Common Name	Count	Cumulative %
Portable XR Chest AP single view	41057	30.7%
XR Chest AP	21416	46.7%
XR Chest Single view	19990	61.6%
XR Chest 2 Views	9754	68.9%
XR Chest PA and Lateral	9173	75.7%
Portable XR Chest Views AP	6533	80.6%
XR Unspecified body region Views	6494	85.4%
CT Chest WO contrast	4148	88.5%
CT Abdomen W contrast IV	115	99%

Row 22

MIDRC-LOINC Mapping Table

DICOM		LOINC	
Modality	Study Description	LOINC code	LOINC Long Common Name
CT	CT CHEST WITHOUT CONTRAST	29252-4	CT Chest WO contrast
CT	CHEST WITHOUT IV CONTRAST	29252-4	CT Chest WO contrast
CT	CT CHEST WO	29252-4	CT Chest WO contrast
CT	CT CHEST WO CONTRAST	29252-4	CT Chest WO contrast
CT	CT CHEST WO IV CONT	29252-4	CT Chest WO contrast
CT	CT CHEST WO IV CONTRAST	29252-4	CT Chest WO contrast
CT	Thorax^CHEST_WITHOUT (Adult)	29252-4	CT Chest WO contrast
CT	Thorax^CHEST_WO (Adult)	29252-4	CT Chest WO contrast
CT	Thorax^CT_CHEST_WITHOUT (Adult)	29252-4	CT Chest WO contrast
CT	Thorax^CHEST_WO_GR (Adult)	29252-4	CT Chest WO contrast
CT	Thorax^ROUTINE_CHEST_WO (Adult)	29252-4	CT Chest WO contrast

The MIDRC-LOINC Mapping Table is open source on Github

Using LOINC Codes for Cohort Building

- Say you find a Study Description “Thorax^CHEST_WO (Adult)” that looks good
- In the MIDRC Data Explorer table, the LOINC Long Common Name is “CT Chest WO contrast”
- Selecting the Long Common Name filter tab, search for “CT Chest WO contrast”
- This brings up over 700 studies, with about 20 different Study Descriptions like:
 - Thorax^XL_CHEST_WO (Adult)
 - CT CHEST WO
 - CT CHEST WO CONTRAST
 - CT CHEST WO IV CONT
 - CT CHEST HIGH RESOLUTION
 - CT CHEST WITHOUT IV CONTRAST

Inspecting Studies

- We are not necessarily done yet
 - Some studies may be mis-labelled by a submitting site
 - Some useful studies may not be in the selected cohort
 - We can inspect images directly

The screenshot displays a web application for managing imaging studies. On the left, a sidebar contains a 'Filters' section with tabs for 'Annotations', 'Case Demographics', 'Imaging Properties', 'LOINC Properties', and 'Study Properties'. The 'Study Properties' tab is active, showing a 'Study UID' dropdown and a 'Study' filter with 'CT' selected. The main area shows a table of 'Imaging Studies' with 735 total entries. The table has columns for 'Study Modality', 'Study Description', and 'Body Part Examined'. The first row shows 'CT' modality, 'Thorax^XL_CHEST_WO (Adult)' description, and 'CHEST' body part. A red box highlights the 'Browse in DICOM viewer' button in the first row, with a red arrow pointing to it and the text 'Click on this'.

Filters

Annotations

Case Demographics

Imaging Properties

LOINC Properties

Study Properties

[Collapse all](#)

Study UID

Select...

Study

1 selected

CT

735

Download Table

Download File Manifest for Imaging Studies (4.86K)

Imaging Studies

735


Showing 1 - 20 of 735 imaging_studies

Show Empty Columns

	Study Modality	Study Description	Body Part Examined
Browse in DICOM viewer	CT	Thorax^XL_CHEST_WO (Adult)	CHEST

Click on this

- This launches the OHIF viewer
- Many useful tools, including a DICOM meta-data viewer
- All the series are listed and can be viewed



Downloading Studies

- Whatever filters are applied selects the cohort manifest that will be downloaded
- The gen3-client application uses the manifest to pull the images locally

Click on this



Filters

Annotations

Case Demographics

Imaging Properties

LOINC Properties

Study Properties

[Collapse all](#)

▼ Study UID

Select... ▼

▼ Stu

1 selected

✕

⚙

🔍

☒ CT

735

Download Table ⋮


Download File Manifest for Imaging Studies (4.86K) ⬇

Imaging Studies

735

Showing 1 - 20 of 735 imaging_studies

Show Empty Columns ☐

Browse in DICOM viewer	Study Modality	Study Description	Body Part Examined
	CT	Thorax^XL_CHEST_WO (Adult)	CHEST

Summary

- LOINC Long Common Names can provide a more meaningful and systematic set of study descriptions
- MIDRC uses a subset (currently < 200) of LOINC codes to simplify searching
- Some details may be lost, as there are almost 2k study descriptions in MIDRC, and over 7k radiologic LOINC codes
- This is due to the level of detail sometimes needed in local clinical use of a study description, versus what is desired to build a cohort from pooled data
- Even with the use of the MIDRC-LOINC descriptions, some data wrangling may be required
- You can use the MIDRC-LOINC descriptions that are built into the MIDRC data explorer, or access directly from Github
- Other LOINC filters: Contrast, Method (i.e. modality), and System (i.e. body region)

Links

- MIDRC data explorer data.midrc.org/explorer
- MIDRC-LOINC mapping table for study descriptions on [Github](#)
- Radlex Playbook / LOINC codes loinc.org (use “class:rad” in LOINC search)